Product datasheet

Quadroflex Advanced

A 4-section couch; as versatile as you. With hinges beneath your patient's knee, hip and shoulder joints, the Quadroflex provides ideal support. This anatomically based split makes almost every position possible. Ideal in case of moving your patient from laying to sitting position. The patient does not have to rotate his/her feet on the head section to take a sitting position which raises comfort and hygiene.

Positioning with the Quadroflex is easy thanks to the 3 powerful motors.

Standard features

- Four-section couch with adjustable head-, back-, sitting-, and leg section
- A total of 3 motors control height, back rest and sitting position
- Arm supports in the head section
- Optional fold-down supports in the back section
- Electric height adjustment, operated by
- peripheral foot switch from all 4 sides of the couch, or single footswitch
- Gymna Stability Profile (GSP)
- Cushions with 40 mm padding
- Breathing hole in head section, incl. cushion
- Central wheel-lifting mechanism with 4 casters
- Non slip feet
- Cushions in a choice of 11 colours
- Prepared for shoulder traction supports

Ordering data

- 125 091 Quadroflex Advanced peripheral switch
- 125 092 Quadroflex Advanced foot switch

Length x width:
Height adjustment:
Weight:
Max. lifting capacity:
Electric motor:
Power supply:
Current consumption:
Insulation class:
Moister resistance:
Classification:
Conformity:

± 202 x 67cm ± 45 – 95 cm ± 110 kg 200 kg 6.000 N 100 – 240 VAC, 50-60 Hz 1,5 A II IP 21 I IEC/EN 60601-1 IEC/EN 60601-2 CAN/CSA C22.2 No. 601.1 MDD 93/42/EEC

3 years (upholstery 1 year)



Guarantee:

Optional

125 618 Fold-down sections (Advanced)

Accessories

- 122 405 Paper roll holder, for undercarriage 337 018 Paper rolls (per 12 pieces) 107 819 Traction support for shoulders 121 401 Gynaecological leg support 333 245 Traction device AutoTrac 460 121 434 Attachment frame for AutoTrac 450/460 332 915 Foot end protector for cushion XXX XXX Terry-cloth cover, stretch material
 - (available in 9 colours 'Style')



Gymna Stability Profile



